

Gang Tensile Core Box



1.0 Product Information:

1.1 General Information

The equipment consists of split core box, base plate and allen screw for clamping arrangement.

1.2 General Usage

This equipment used for prepare the tensile type multiple specimen of No-bake sand.

1.3 Product Data

Product Name:	Gang Tensile Core Box
Product Code:	VGT
Order Reference Code:	V0000020
Electrical Supply Required:	Not Required
Electrical Plug:	Not Applicable
Typical Maximum Current:	Not applicable
Compressed air Supply:	Not Required
Communication:	Not Applicable
Dimensions:	260x185x25mm (WXBXH)
Packing Dimensions:	350x250x150mm (WXBXH) Typical
Net Weight:	Approximately 3 kg
Gross Weight:	Approximately 4 kg
Colour:	Factory Default
Measuring Range:	Not Applicable
Traceable to:	Manufacturers Standards
Least Count:	Not Applicable
Specimen Requirements:	Not Required
Operating Temperature:	0-50°C
Testing Standards:	BIS: Standard AFS: Optional GOST: Optional
Calibration Frequency:	4500 nos. of Testing Cycles or Once in Year

2.0 Scope of Supply

2.1 Inclusions in Supply Scope

- Tensile core box.
- Base plate with allen screws.
- Instruction manual.

2.2 Exclusions in Supply Scope

- Anything that is not covered under Supply scope.

3.0 Recommended Pre-requisite Equipment

NA

4.0 Accessories

4.1 Calibration Kits

NA

4.2 Attachments

NA

5.0 Alternative Equipment

[Tensile Core Box VCB](#) ↗

6.0 Installations Instructions

Required Levelled (2mm) and sturdy platform. Other Installation instruction shall be provided along with the equipment.

7.0 Consumables

NA

8.0 Hazards

No Major Hazards.

9.0 Additional Information

[Link to Web-Page](#) ↗

Disclaimer: Specifications and design of equipment may change without prior notice.